




# NEW HORIZONS — ALLERGY —

*New recombinant markers for use in immunotherapy*  
Join us exploring the frontiers of allergy research



The use of these new recombinant allergen markers in immunotherapy may:

- assist in inclusion of suitable patients (IgE)
- demonstrate induction of immunological response (IgG)
- reveal the development of cross-reactive IgE specificities

**rPar j 2** represents a marker allergen for genuine sensitization to *Parietaria* pollen and can thus be used for the same concept.



## New Horizon Products

- Rg213 rPhl p 1, rPhl p 5b
- Rg214 rPhl p 7, rPhl p 12
- Rw211 rPar j 2
- Rt215 rBet v 1
- Rt221 rBet v 2, rBet v 4

New recombinant markers for use in immunotherapy.  
Join us exploring the frontiers of allergy research!



## Product List

### ImmunoCAP™ IgE inclusion tests for immunotherapy

#### Grass pollens

Code	Product	Art. No.
Rg213	rPhl p 1, rPhl p 5b	14-5312-10
Rg214	rPhl p 7, rPhl p 12	14-5313-10

#### Birch pollens

Code	Product	Art. No.
Rt215	rBet v 1	14-5225-10
Rt221	rBet v 2, rBet v 4	14-5310-10

#### Parietaria pollens

Code	Product	Art. No.
Rw211	rPar j 2	14-5311-10

### ImmunoCAP™ IgG follow-up tests for immunological response

#### Grass pollens

Code	Product	Art. No.
Rg213	rPhl p 1, rPhl p 5b	14-5312-10

#### Birch pollens

Code	Product	Art. No.
Rt215	rBet v 1	14-5225-10

#### Parietaria pollens

Code	Product	Art. No.
Rw211	rPar j 2	14-5311-10

### ImmunoCAP™ IgE follow-up tests for cross-reactive specificities

#### Grass pollens

Code	Product	Art. No.
Rg214	rPhl p 7, rPhl p 12	14-5313-10

#### Birch pollens

Code	Product	Art. No.
Rt221	rBet v 2, rBet v 4	14-5310-10

**Phadia**

Phadia AB, P O Box 6460, SE-751 37 Uppsala, Sweden  
Tel +46 18 16 50 00. [www.phadia.com](http://www.phadia.com)